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Mina Bone Master All-around Kit is a multi-functional dental bone tissue kit suitable for performing operations such as bone extraction, grinding, and cutting in dental surgeries. It includes various types of drills, such as Pillar trephine drill, precision drill (for positioning), round bur, bone trimmer, safe click and bone shaver etc, which can meet different surgical needs.

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1. According to the needs of the surgery, select appropriate size and diameter of the pillar trephine drill.
2. Use the pillar trephine drill to extract a columnar bone block in the implant site.
3. Transplant the extracted columnar bone block to bone defect position.
4. For small bone defects around the implant site, it can be used directly without additional use of bone powder or membranes.
5. The removed columnar bone block can also be crushed and mixed with bone powder for use, usually recommended at a 1:1 ratio.
6. During the use, by using the special detachable working end and texture design of the pillar trephine drill, it is convenient to extract the entire columnar bone.
7. Through the built-in internal cooling port of the drill, double flushing is adopted to protect the bone tissue from overheating and necrosis.

产品构成Components

01 取骨柱钻 (Pillar Trepine Drill)

- 取骨柱钻壁厚0.5mm，内径-1mm
- 取骨柱钻长38mm(HL)
- 可取骨柱长度16mm(L)
- 推荐转速:600~1000 rpm
- The wall thickness of Pillar Trepine Drill is 0.5mm, and the outer diameter is 1 mm larger than inner diameter.
- Total length of the Pillar Trepine Drill is 38mm(HL).
- The bone column can reach 16mm(L).
- Recommended speed: 600~1000rpm.

型号 Model	BPT21	BPT25	BPT29	BPT33
				
D (mm)	3.1	3.5	3.9	4.3
L (mm)	16	16	16	16
HL (mm)	38	38	38	38

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It's applicable for extracting the bone columns with required diameter.

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
The pillar trephine drill is detachable for easy removal of bone column.


1. Use the HRK hand wheel to clamp the bottom(implanter connection) and rotate the drill part.
2. After detached, put the drill into the ASH ejection sleeve from top.
3. Push out the bone column with applicable ejection rod.
4. Stuck bone can be removed with AS08 auxiliary rod.

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02 定位取骨钻 (Guide Bone Trephine)

- 定位取骨钻壁厚0.5mm，内径-1mm
- 定位取骨钻长22mm(HL)
- 定位深度4.2mm(L)
- 推荐转速:800~1000rpm
- The wall thickness of guide bone trephine is 0.5mm, and the outer diameter is 1 mm larger than inner.
- Total length of guide bone trephine is 22mm(HL).
- Guide depth is 4.2mm(L).
- Recommended speed: 800~1000rpm.


型号 Model	GBT21	GBT25	GBT29	GBT33
				
D (mm)	3.1	3.5	3.9	4.3
L (mm)	4.2	4.2	4.2	4.2
HL (mm)	22	22	22	22

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It's used for the preliminary positioning before the Pillar Trephine Drill is applied.


03 环钻辅助套(Ejection Sleeve)

- 环钻辅助套直径12mm,内环直径5.5mm(D)
- 环钻辅助套长28mm(HL)
- 用于配合取骨柱钻、骨柱推动杆、环钻辅助杆取出内部骨柱
- Diameter of ejection sleeve is 12mm, and inner ring is 5.5mm(D).
- Total length of ejection sleeve is 28mm(HL).
- To cooperate with pillar trephine drill, ejection rod and auxiliary rod to remove the internal bone column.

04 骨柱推动杆(Ejection Rod)


- 骨柱推动杆工作端直径分别2.1mm(D)和2.9mm(D)
- 骨柱推动杆工作端长度22mm(L)
- 骨柱推动杆长45.4mm(HL)
- There are 2.1mm and 2.9mm(D) diameter for working end of ejection rods.
- Length of ejection rod working end is 22mm(L).
- Total length of ejection rod is 45.4mm(HL).

型号 Model	AS21	AS29
D (mm)	Ø 2.1	Ø 2.9
L (mm)	22	22
HL (mm)	45.4	45.4

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


It is applicable for preliminary positioning and drilling, and can also be used to drill nutrient foramina to provide blood supply.

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It is applicable for precisely grinding bone tissue and cleaning granulation tissue.

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- The diameter of bone collector BC50 is 5 mm(D).
- Total length of bone collector is 35 mm(HL).
- Length of stopper is 15 mm(L2).
- Recommended speed: 800-1000rpm.

型号 Model	BC50
D (mm)	Ø 5.0
L (mm)	21
HL (mm)	35

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
It's suitable for obtaining autogenous bone chips and has a self-grinding function.
It's used in combination with stopper.

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- The diameter of bone trimmer is 6 mm, and minimum diameter is 4mm(D).
- Working end length is 12mm(L).
- Total length is 29.6 mm(HL).
- Recommended speed: 800-1000rpm.

型号 Model	PN60
D (mm)	Ø 6.0
L (mm)	12
HL (mm)	29.6

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The head is smooth, which can avoid damaging soft tissues during grinding.

The drill body is designed with a concave surface, which can directly trim the alveolar ridge to a regular shape.

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
It's suitable for grinding irregular alveolar ridges.

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

It's applicable for incising the bone surface. It comes with a protective sleeve to avoid damaging the surrounding soft tissues.

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It is suitable for collecting autogenous bone chips. With a detachable blade design, it is convenient for replacement.

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

- ① After using all the tools, immerse them in the disinfectant solution immediately for a few minutes to avoid chips sticking to the tools.
Note: Do not soak them overnight. Do not use hydrogen peroxide for soaking, as hydrogen peroxide will cause the tool markings to change color.
- ② Please use a soft-bristled brush to scrub the tools carefully until all chips are removed.
- ③ Please use a drill or a measuring injection to clean the inside of drill bits.
- ④ If an ultrasonic cleaner is used, please wrap the drill bits with clean gauze before cleaning to avoid friction between the tool drill bits.
- ⑤ Please use disinfectant solution to clean all bases and lids.
- ⑥ Please thoroughly and carefully clean all tools, trays and lids in warm water.
- ⑦ After complete drying, please check carefully to confirm that there is no damage or corrosion on the tools.
- ⑧ Please put all the tools back onto the trays in the toolbox according to position markings.
- ⑨ Please seal the toolbox with a sterilization bag.
- ⑩ Please sterilize the toolbox in a steam autoclave at 134 °C for 15 minutes, or refer to other suggestions from the manufacturer.
- ⑪ Please store the toolbox in a dry and clean area at room temperature.

Precautions:

After the operation, please separate and clean the used instruments in time. It is recommended to disinfect the surgical instruments once again before next operation (temperature: 134 °C, time: 15 minutes).

Company Name: Mapleisland Medical Technology (Qingdao) Co., Ltd.
Address: 77 Keyun Rd, Hi-tech Zone, Qingdao City, Shandong Prov., China
Postal Code: 266114
Mobile: +86-15308385888
Office Tel: +86-532-80963040
Email : mapleisland2024@163.com
WS: www.mnmina.com

Made in China