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
MINA Guide Collector Kit

Edition 2025

定位导向钻/Guide Drill

型号/Model	直径/Diameter (D)	颜色/Color
GD 06	Ø 6.0	水晶紫Crystal violent
GD 07	Ø 7.0	摩登金Modern gold
GD 08	Ø 8.0	宝石绿Jewelry green
GD 09	Ø 9.0	天际蓝Sky blue
GD 10	Ø 10.0	极光蓝Aurora blue

特点Characteristics:

- 1.在备洞的同时,可以收集到自体骨,用于自体骨移植
 - 2.便于在两颗种植体同时种植时,确认两颗种植体的间距,减少手术中的失误
 - 3.简单有效的把握如何确定种植窝位置
- 1.1.Collect autologous bone for autologous bone transplantation while preparing hole.
 - 2.Convenient to confirm distance between two implants when simultaneously implanting two implants, reducing errors during operation.
 - 3.Simple and effective to grasp how to determine the implant socket position.

定位杆/Guide Pin

型号/Model	直径/Diameter(D)	颜色/Color
GP 06	Ø 6.0	水晶紫Crystal violent
GP 07	Ø 7.0	摩登金Modern gold
GP 08	Ø 8.0	宝石绿Jewelry green
GP 09	Ø 9.0	天际蓝Sky blue
GP 10	Ø 10.0	极光蓝Aurora blue

特点Characteristics:

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1. 1.Fixing a single guide pin can confirm the implant socket position.
- 2.It makes surgery more convenient when implanting two or more implants simultaneously.

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
1. 1.After horizontally shaving the ridge top of insufficient width, the position of implant socket can be marked on trimmed alveolar bone plane.

2.When used after extraction of residual tooth, it can effectively remove remaining bone and facilitate formation of implant position.

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1.The autologous bone collected during cavity preparation can be used in bone grafting, thereby avoiding additional autologous bone collection surgeries or using artificial bone powder.
2.Reduce surgical steps and shorten time.

操作实例Operation examples:


▲开始取骨时应将停止环下缘与
钻尖齐平, 可以获得大量自体骨。
The lower edge of stopper should be aligned with
drill tip when starting to collect a large amount of
autogenous bone.

▲开始取骨
Start collecting

▲取骨结束
End of collecting

▲取下停止环, 单次取骨即可
获得大量自体骨
Take off the stopper. It has collected a large
amount of autogenous bone for single time.

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
△1个新型取骨钻
△1 bone collector

规格 Specification (Ø6/7/8/9/10mm)

*所有钻头的直径都是2mm
The diameters of all the drills are 2mm.

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After use, separate guide drill as left picture, and clean it with ultrasonic cleaner.
Easy to separate by pinching it with three fingers.

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If the drill is rotated after contacting alveolar bone, it is prone to slide. The optimum rotational speed (1000-1200rpm) should be reached first before contacting alveolar bone.

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

For preliminary cavity preparation and collecting autologous bone, a rotational speed of 1000-1200rpm should be selected and water injection be performed simultaneously. Water injection can prevent bone heating, but the collected bone may decrease.

直径选择 (Diameter selection)


▲ 如果选择了不合适的定位钻(8mm)会被邻牙卡住
If choose inappropriate drill (8mm), it will be blocked by adjacent tooth

▲ 使用适当的定位钻(7mm)钻尖接触到牙槽骨
Drill tip can touch alveolar bone by using appropriate drill (7mm)

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

Hold the handpiece with right hand and use left hand to stabilize the head of the handpiece when drilling, which can prevent positional deviation caused by handpiece vibration.

*Those who are accustomed to using left hand should use it in the opposite way.

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
Use no. 8 guide drill (Jewelry green) to determine #45 position. The max drill length is 7.2mm, which can be referred to determine drilling depth. The guide drill can collect autologous bone while drilling.

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
Insert no. 8 guide pin at the #45 position and the 8# guide pin represents the second premolar now. Use no. 9 guide drill to accurately locate the #46 position.

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
Insert no. 8 and 9 guide pin at #45 and #46 to mark the second premolar and first molar.

To determine correct position of #47, a no. 10 guide drill should be used to drill, and the drill should be parallel to adjacent guide pins.

Good contact with adjacent guide pins is important for guide.

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

As above picture, #45, #46 and #47 bite the guide pin corresponding to the maxillary teeth, which can determine the occlusal point and facilitate accurate repair in later stage.

#45-#47种植牙手术 (Dental Implant Surgery)

钻孔完成后, 植入种植体。植入深度参考定位杆的“使用方法”匹配使用。
The implant is implanted after drilling. The depth refer to “Usage Method” of the guide pin for matching use.

定位导向钻
Guide Drill
($\varnothing 6$ 、 $\varnothing 7$ 、 $\varnothing 8$ 、 $\varnothing 9$ 、 $\varnothing 10$)

定位取骨工具箱
Guide Collector Kit

定位杆
Guide Pin
($\varnothing 6$ 、 $\varnothing 7$ 、 $\varnothing 8$ 、 $\varnothing 9$ 、 $\varnothing 10$)

$\varnothing 5.0$ 新型取骨钻
 $\varnothing 5.0$ Bone Collector

$\varnothing 5.0$ 菠萝钻
 $\varnothing 5.0$ Bone Trimmer

Edition 2025

Company Name: Shandong YuFeng JiaYe Medical Technology Co., Ltd.
Address: 77 Keyun Rd, Hi-tech Zone, Qingdao City, Shandong Prov., China
Postal Code: 266114
Mobile: +86-15308385888
Office Tel: +86-532-80963040
Email : mapleisland2024@163.com
WS: www.mnmina.com